




Standard Extension Set 80cm

PRODUCT DATA SHEET



DESCRIPTION

The m|devices Standard Extension Set is 80cm in length with two NeutralSite™, two slide clamps and a female to male luer lock line with rotating collar (RC).


FEATURES AND BENEFITS

- 2 x NeutralSite™
- 2 x Slide clamps
- Standard bore tubing
- Transparent, soft, smooth tubing
- Female luer lock to male luer lock end with RC
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 80cm
- Priming volume: 6.01mL
- NeutralSite™: Neutral displacement swabable needleless access site
- Sterile EO
- MRI compatible
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Not made with Natural Rubber Latex
- Material: HDPE, PC, PE, DEHP Free PVC tubing, ABS, LDPE
- NeutralSite™ Materials: PP, ABS and Silicone.
- Shelf Life: 3 Years
- ARTG Number: 354371

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV052201	Standard Extension Set	80cm	100

m|devices®

m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au