




# Standard Extension Set

## 80cm

### PRODUCT DATA SHEET



#### DESCRIPTION

The m|devices Standard Extension Set is 80cm in length with two swabable Needleless Access Sites, two slide clamps and a female to male luer lock line with rotating collar (RC).


#### FEATURES AND BENEFITS

- 2x Needleless Access Site
- 2 Slide clamps
- Standard bore tubing
- Transparent, soft, smooth tubing
- Female luer lock to male luer lock end with RC
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only

#### PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 80cm
- Priming volume: 6.01mL
- DEHP free
- Sterile EO
- Latex free
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- MRI compatible

#### PRODUCT PACKAGING LABEL



| RE-ORDER | DESCRIPTION            | SIZE | CARTON/QTY |
|----------|------------------------|------|------------|
| IV052201 | Standard Extension Set | 80cm | 100        |



**m|devices**  
27 Llewellyn Avenue Villawood  
NSW 2163 Australia  
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

[info@mdevices.com.au](mailto:info@mdevices.com.au)  
[mdevices.com.au](http://mdevices.com.au)