

Standard Extension Set 80cm

DESCRIPTION

The m|devices Standard Extension Set is 80cm in length with two NeutralSite™, two slide clamps and a female to male luer lock line with rotating collar (RC).

FEATURES & BENEFITS

- 2 x NeutralSite™
- 2 x Slide clamps
- Standard bore tubing
- Transparent, soft, smooth tubing
- Single use only
- Female luer lock to male luer lock end with RC
- Made from non-toxic medical grade PVC
- Peel open pouch


PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 80cm
- Priming volume: 6.01mL
- Sterile EO
- NeutralSite™: Neutral displacement swabable needleless access site
- ARTG Number: 354371
- MRI compatible
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Material: HDPE, PC, PE, DEHP Free PVC tubing, ABS, LDPE
- NeutralSite™ Materials: PP, ABS and Silicone
- Shelf Life: 3 Years

ISO STANDARDS

- ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment
—Part 1: General requirements
- ISO 594-2 Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment
—Part 2:Lock fittings
- ISO14971 Medical Devices- Application of risk management to medical device
- ISO13485 Medical Devices- Quality management systems – Requirements for regulatory purposes

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Clamp
4	Tubing
5	NeutralSite™
6	Male Luer Lock
7	Cap

Standard Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052201	Standard Extension Set	20cm	100