

Standard Extension Set

120cm

DESCRIPTION

The m|devices Standard Extension Set is 120cm in length with two Luer-Site needleless sites, two slide clamps and a female to male luer lock line with rotating collar (RC).

FEATURES & BENEFITS

- 2 x Luer-Site needleless access site
- 2 x Slide clamps
- Standard bore tubing
- Transparent, soft, smooth tubing
- Female luer lock to male luer lock end with RC
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: I.D- 3mm O.D- 4mm
- Tubing length: 120cm
- Priming volume: 8.48mL
- Sterile EO
- MRI compatible
- Lipid resistant: No
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 354371

ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Female Luer Lock
2	Slide Clamp
3	Luer-Site Needleless Access Site
4	Tubing
5	Male Luer Lock

Standard Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052203	Standard Extension Set	120cm	100