




Standard Extension Set

120cm

PRODUCT DATA SHEET



DESCRIPTION

The m|devices Standard Extension Set is 120cm in length with two Luer-Site needleless sites, two slide clamps and a female to male luer lock line with rotating collar (RC).


FEATURES AND BENEFITS

- 2 x Luer-Site needleless access site
- 2 x Slide clamps
- Standard bore tubing
- Transparent, soft, smooth tubing
- Female luer lock to male luer lock end with RC
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 120cm
- Priming volume: 8.48mL
- DEHP free
- Sterile EO
- Lipid resistant: No
- Latex free
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- MRI compatible

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV052203	Standard Extension Set	120cm	100



m|devices
27 Llewellyn Avenue Villawood
NSW 2163 Australia
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
mdevices.com.au