




# Multi Lumen Standard Bore Extension Set

15cm  
2 Lumens

## PRODUCT DATA SHEET



### DESCRIPTION

The m|devices Multi Lumen Standard Bore Extension Set with two Lumens, is 15cm in length. It has a female luer lock and slide clamp on each lumen and a male luer lock connector and rotating collar at the other end.


### FEATURES AND BENEFITS

- 2 x Female luer lock connectors
- 2 x slide clamps
- 2 x lumen
- Male luer lock with rotating collar
- Peel open pouch
- Single use only

### PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 15cm
- Priming volume (2 lumen): 1.60mL
- DEHP free
- Sterile EO
- Latex free
- Material: PVC, PE, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- The product is MRI and CT suitable

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV052300	Multi Lumen Standard Bore Extension Set	15cm	100

m|devices®

m|devices  
27 Llewellyn Avenue Villawood  
NSW 2163 Australia  
T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au  
mdevices.com.au