

Multi Lumen Standard Bore Extension Set

15cm
3 Lumens

DESCRIPTION

The Multi Lumen Standard Bore Extension Set with three Lumens, is 15cm in length. It has female luer lock and slide clamp on each lumen and a male luer lock connector and rotating collar at the other end.

FEATURES & BENEFITS

- 3 x Female luer lock connectors
- 3 x slide clamps
- 3 x lumens
- Male luer lock with rotating collar
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: I.D- 3mm O.D- 4mm
- Tubing length: 15cm
- Priming volume (3 lumen): 4.10mL
- Lipid resistant
- Sterile EO
- ARTG Number:354371
- The product is MRI and CT suitable
- Not made with DEHP
- Not made with Natural Rubber Latex
- Material: HDPE & LDPE Caps, PE side clamp, PVC Tubing, PVC connector, ABS MLL.
- Shelf Life: 3 Years

ISO STANDARDS

- ISO 13485-2016 Quality management systems. Requirements for regulatory purposes
- ISO 1497-2012 Medical devices. Application of risk management to medical devices
- ISO 11135-2014 Sterilisation of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilisation process for medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Female Luer Lock
3	Slide Clamp
4	Tube
5	Connector
6	Tube
7	Male Luer Lock
8	Cap

Multi Lumen Standard Bore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052400	Multi Lumen Standard Bore Extension Set	15cm	100