

Standard Extension Set 15cm

DESCRIPTION

The m|devices Standard Extension Set is a 15cm tubing with a negative pressure needleless access site through to a male luer lock connector with rotating collar (RC).

FEATURES & BENEFITS

- Standard Bore tubing
- Needleless access site
- Male luer lock end with rotating collar
- Transparent, soft, smooth tubing
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only


PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 15cm
- Priming volume: 1.3mL
- Sterile EO
- Material: PE, ABS, PVC, PC Silicone
- Not made with DEHP
- Not made with Natural Rubber Latex
- Shelf Life: 3 Years
- ARTG Number: 354371

ISO STANDARDS

- ISO 8536-4: Infusion equipment for medical use- Part 4: Infusion sets for single use, gravity feed
- ISO 80369-7: Small-bore connectors for liquids and gases in healthcare applications-Part 7: Connectors for intravascular or hypodermic applications
- ISO14971: Medical Devices- Application of risk management to medical device

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Male Luer Lock
3	Tubing
4	Female Luer Lock
5	Needleless Access Site

Standard Bore Extension Set

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV052416	Standard Bore Extension Set	15cm	200