




# Microbore Extension Set

## 15cm

### PRODUCT DATA SHEET



#### DESCRIPTION

The Microbore Extension Set is 15cm in length with two lumens, each with a swabable Needleless Access Site and slide clamp to a distal male luer lock connector with a rotating collar (RC).


#### FEATURES AND BENEFITS

- 2 x Lumens, each with a Needleless Access Site
- Microbore tubing
- 2 x Lumen slide clamps
- Male luer lock with rotating collar
- Transparent, soft, smooth tubing
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only

#### PRODUCT DATA

- Tubing diameter: 1mm
- Tubing length: 15cm
- Priming volume: 0.17mL
- DEHP free
- Sterile EO
- Lipid resistant: No
- Latex free
- Material: PVC, PE, PC
- Shelf Life: 5 Years
- ARTG Number: 228230
- MRI compatible

#### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV053001	Microbore Extension Set	15cm	100



**m|devices**  
 27 Llewellyn Avenue Villawood  
 NSW 2163 Australia  
 T + 61 (0)2 8718 2800 F + 61 (0)2 8718 2801

[info@mdevices.com.au](mailto:info@mdevices.com.au)  
[mdevices.com.au](http://mdevices.com.au)