



## Minimum Volume Extension Set 15cm

### PRODUCT DATA SHEET



#### DESCRIPTION

The m|devices Minimum Volume Extension Set is a 15cm tubing with a swabable Needleless Access Site and lumen click clamp through to a male luer lock connection end with rotating collar (RC).



#### FEATURES AND BENEFITS

- Needleless Access Site
- Microbore tubing
- Lumen click clamp
- Male luer lock end with RC
- Transparent, soft, smooth tubing
- Made from non-toxic medical grade PVC
- Peel open pouch
- Single use only

#### PRODUCT DATA

- Tubing diameter: 1.5mm
- Tubing length: 15cm
- Priming volume: 0.26mL
- Sterile EO
- Lipid resistant: No
- PSI rating: 325psi
- DEHP free
- Latex free
- Material: PVC, PE, ABS, PC
- Shelf Life: 3 Years
- ARTG Number: 228230
- MRI compatible

#### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QT
IV054000	Minimum Volume Extension Set	15cm	200