

IV Transfer Set

Vented 20cm

PRODUCT DATA SHEET



DESCRIPTION

The IV Transfer Set has a vented perforator and click clamp at one end and a male luer lock at the other.

FEATURES AND BENEFITS

- Vented perforator
- Male luer lock with rotating collar
- Click clamp
- Peel open pouch
- Sterile EO
- Latex free
- DEHP free

PRODUCT DATA

- Tubing length: 20cm
- Tubing diameter: 3mm
- Priming volume: 2.0mL
- Material: PVC
- Shelf Life: 5 Years
- ARTG Number: 160343
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
IV070000	20cm	100

RE-ORDER	SIZE	CARTON QTY