



Double Spike Adaptor 30cm

PRODUCT DATA SHEET



DESCRIPTION

The Double Spike Adaptor is designed with two vented perforators with roller regulators and a sealed septum connector. Made from medical grade non-toxic PVC (latex free).

FEATURES AND BENEFITS

- 2 x Vented perforators
- 2 x Roller regulators
- Sealed septum connector
- Allows adaption from single to bag to double bag
- Peel open pouch

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 30cm
- Priming volume: 5.4mL
- Sterile EO
- DEHP / Latex free
- Material: PVC
- Shelf Life: 5 Years
- ARTG Number: 160343
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV071000	Double Spike Adaptor	30cm	100