

Double Spike Adaptor 30cm

DESCRIPTION

The Double Spike Adaptor is designed with two perforators with roller regulators and a septum connector. Made from medical grade non-toxic PVC (latex free).

FEATURES & BENEFITS

- 2 x Perforators
- 2 x Roller regulators
- Septum connector
- Allows adaption from single to bag to double bag
- Peel open pouch


PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 30cm
- Priming volume: 5.4mL
- Sterile EO
- The product is MRI and CT suitable
- ARTG Number: 304860
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected [Result less than the Method Detection Limit EN 14372]
- Material: PVC
- Shelf Life: 3 Years

ISO STANDARDS

- ENISO8536-4:2013 Infusion equipment for medical use - part4: infusion set
- ISO14971:2016 Medical Devices - Application of risk management to medical device
- ISO13485:2016 Medical Devices - Quality management systems- Requirements for regulatory purposes

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Spike
3	Tubing
4	Roller Clamp
5	Connectors
6	Tubing
7	Connectors

Double Spike Adaptor

STERILE

RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
IV071000	Infusion/Transfusion Double Spike Adaptor	30cm	50