




Closed System Transfer Set

43cm

PRODUCT DATA SHEET



DESCRIPTION

The Closed System Transfer Set is a needle free device with multi-lumen extensions designed for efficient fluid transfer. It consists of 2x2 lumen female luer lock connectors, attached to a three-way stopcock. It also features three perforators with click clamps and a swabable needleless access valve.


FEATURES AND BENEFITS

- Dedicated set for efficient fluid transfer
- 2 x 2 Lumen female luer lock connectors
- 3 - Way stopcock
- 3 x Vented perforators
- 1 x swabable Needleless Access Site
- Long length with multiple points
- Peel open pouch

PRODUCT DATA

- Tubing diameter: 3mm
- Tubing length: 43cm
- Priming volume: 5mL
- Temperature range: 4 - 60°C
- The product is MRI and CT suitable
- Sterile EO
- DEHP free
- Latex free
- Material: PVC, PC, ABS, PE
- Shelf Life: 3 Years
- ARTG Number: 228230

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV072001	Closed System Transfer Set	43cm	20



m|devices
 27 Llewellyn Avenue Villawood
 NSW 2163 Australia
 T + 61 (0) 2 8718 2800 F + 61 (0) 2 8718 2801

info@mdevices.com.au
 mdevices.com.au