

3 Way Stop Cock

PRODUCT DATA SHEET



DESCRIPTION

The 3-Way Stop Cock allows for 360° rotation with a transparent housing for greater visibility. It requires only a minimum amount of priming fluid and can take pressures up to 4.48bar (65Psi). The arrows on the handle indicate the direction of flow.

FEATURES AND BENEFITS

- 360° rotation
- Transparent channel
- Minimum priming volume required
- Male luer lock with rotating collar
- Two female luer connectors
- Pressure up to 4.48bar (65Psi)
- Arrow on handle indicates direction of flow
- DEHP free

PRODUCT DATA

- Priming volume: 0.3mL
- Shelf life: 5 Years
- Material: Polycarbonate
- ARTG Number: 278076
- Colour: White
- Latex free
- PVC free
- Single use only
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY
----------	------	------------

IV110000		50
----------	--	----