

3-Way Stopcock

Blue
Lipid Resistant



PRODUCT DATA SHEET

DESCRIPTION

The 3-Way Stopcock Lipid Resistant allows for 360° rotation with a transparent housing for greater visibility. It requires only a minimum amount of priming fluid and can take pressures up to 4.48bar (65Psi). The arrows on the handle indicate the direction of flow.


FEATURES AND BENEFITS

- 360° rotation
- Transparent channel
- Minimum priming volume required
- Male luer lock with rotating collar
- Two female luer connectors
- Pressure up to 4.48bar (65Psi)
- Arrow on handle indicates the direction of flow
- Polycarbonate Lipid Resistant material
- Single use only

PRODUCT DATA

- Priming volume: 0.3mL
- Colour: Blue
- DEHP free
- Latex free
- PVC free
- Material: Polycarbonate
- Shelf Life: 5 Years
- ARTG Number: 278076
- The product is MRI and CT suitable

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	COLOUR	CARTON/QTY
IV110001	3-Way Stopcock	Blue	50