




# 3-Way Stopcock with Extension

## PRODUCT DATA SHEET



### DESCRIPTION

The m|devices 3-Way Stopcock with Extension sets come in a variety of lengths; 15cm, 25cm, 50cm, 80cm and 100cm to meet all clinical requirements. Featuring a lipid resistant 3-Way Stop Cock connector to medical grade, non-toxic PVC tubing and make luer lock connector.

### FEATURES AND BENEFITS

- Flow direction arrow indicated on handle
- 3-Way stop cock
- Lipid resistant and blood compatible
- 360° rotation
- Luer slip and luer lock compatible
- Rotating male luer and two female luer connectors
- Peel open pouch

### PRODUCT DATA

- Tubing diameter: 3.0mm
- Priming volume:
 

15cm	1.4mL
25cm	2.1mL
50cm	3.8mL
80cm	5.9mL
100cm	7.3mL
- Flow rate: Up to 174.9mL (gravity)
- Pressure: Up to 4.48bar (45Psi)
- ARTG Number: 228230
- Shelf life: 5 years
- Material: PVC, PC and PP
- DEHP free

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV110002	3-Way Stopcock with Extension	15cm	50
IV110003	3-Way Stopcock with Extension	25cm	50
IV110004	3-Way Stopcock with Extension	50cm	50
IV110005	3-Way Stopcock with Extension	80cm	50
IV110006	3-Way Stopcock with Extension	100cm	50



**m|devices**  
 27 Llewellyn Avenue Villawood  
 NSW 2163 Australia  
 T + 61 (0)2 8718 2800 F + 61 (0)2 8718 2801

[info@mdevices.com.au](mailto:info@mdevices.com.au)  
[mdevices.com.au](http://mdevices.com.au)