

3-Way Stopcock with Extension

DESCRIPTION

The m|devices 3-Way Stopcock with Extension sets come in a variety of lengths; 15cm, 25cm, 50cm, 80cm and 100cm to meet all clinical requirements. Featuring a lipid resistant 3-Way Stopcock connector to medical grade, non-toxic PVC tubing and male luer lock connector.

FEATURES & BENEFITS

- Flow direction arrow indicated on handle
- 3-Way Stopcock
- Lipid resistant and blood compatible
- 360° rotation
- Luer slip and luer lock compatible
- Rotating male luer and two female luer connectors
- Peel open pouch


PRODUCT DATA

- Tubing diameter: ID 3.0mm OD 4.0mm
- Priming volume: 15cm 1.4mL
25cm 2.1mL
50cm 3.7mL
80cm 5.8mL
100cm 7.2mL
- Max pressure: 10 psi
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Materials: DEHP Free PVC, PE, PC, ABS
- ARTG: 354371
- Shelf Life: 3 Years

ISO STANDARDS

- EN ISO 8536-4 Infusion equipment for medical use- Part 4:Infusion sets for single use, gravity feed
- EN ISO 594-1 Conical fittings with a 6% (Luer) taper for syringes, needles and certain other medical equipment —Part 1: General requirements
- EN ISO 594-2 Conical fittings with 6 % (Luer) taper for syringes, needles and certain other medical equipment —Part 2:Lock fittings

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	3-Way Stopcock
2	Tubing
3	Male Luer Lock

3-Way Stopcock with Extension

3-Way Stopcock with Extension			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV110002	3-Way Stopcock with Extension	15cm	50
IV110003	3-Way Stopcock with Extension	25cm	50
IV110004	3-Way Stopcock with Extension	50cm	50
IV110005	3-Way Stopcock with Extension	80cm	50
IV110006	3-Way Stopcock with Extension	100cm	50