

Insulin Syringes

Fixed and Safety

DESCRIPTION

m|devices disposable 1mL syringes have pre-attached needles available in 8mm and 13mm lengths and 27G, 29G and 30G options.

FEATURES & BENEFITS

- Bold scale markings
- Orange cap colour for easy identification
- Pre-attached needle
- Luer slip and luer lock compatible
- Peel open pouch


PRODUCT DATA

- Volume: 1mL and 0.5mL options
- Needle size: 27G, 29G and 30G
- Needle length: 8mm and 13mm
- Needle type: 60-degree bevel
- ARTG: 272214
- Not made with DEHP
- Materials: PP and silicone syringe and Stainless steel needle
- Shelf Life: 5 Years

ISO STANDARDS

- ISO 7886-1:2017 Sterile hypodermic syringes for single use
- EN ISO13485 :2016-Medical devices - Quality management systems - Requirements for regulatory purposes

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Cap
2	Needle
3	Barrel
4	Plunger

Insulin Syringes

Fixed and Safety

Insulin Syringe Fixed				STERILE
RE-ORDER	DESCRIPTION	SIZE	VOLUME	CARTON/QTY
IV120127	Insulin Syringe with Standard Fixed Needle	27G x 13mm	1mL	100
IV120129	Insulin Syringe with Standard Fixed Needle	29G x 13mm	1mL	100
IV120130	Insulin Syringe with Standard Fixed Needle	30G x 13mm	1mL	100
IV120529	Insulin Syringe with Standard Fixed Needle	29G x 13mm	0.5mL	100
IV121127	Insulin Syringe with Standard Fixed Needle	27G x 8mm	1mL	100
IV121129	Insulin Syringe with Standard Fixed Needle	29G x 8mm	1mL	100
IV121130	Insulin Syringe with Standard Fixed Needle	30G x 8mm	1mL	100