




Insulin Syringes

Fixed and Safety
Retractable options

PRODUCT DATA SHEET



DESCRIPTION

m|devices disposable 1mL syringes have pre-attached needles available in 8mm and 13mm lengths and 27G, 29G and 30G options. The safety retractable version retracts the needle into the barrel after injection to prevent re-use and needle stick injury.

FEATURES AND BENEFITS

- Bold scale markings
- Orange cap colour for easy identification
- Pre-attached needle
- Retractable version permanently disables needle reduces risk of needle stick injury and blood splatter
- Luer slip and luer lock compatible
- Peel open pouch

PRODUCT DATA

- Volume: 1mL and 0.5mL options
- Needle size: 27G, 29G and 30G
- Needle length: 8mm and 13mm
- Needle type: 60-degree bevel
- Material: PP and silicone syringe and Stainless-steel needle
- Sterile EO
- DEHP free
- ARTG Number: 272214
- Shelf life: 5 years

RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV120127	Insulin Syringe with Fixed Needle	1mL 27G x 13mm	100
IV120129	Insulin Syringe with Fixed Needle	1mL 29G x 13mm	100
IV120130	Insulin Syringe with Fixed Needle	1mL 30G x 13mm	100
IV120529	Insulin Syringe with Fixed Needle	0.5mL 29G x 13mm	100
IV120530	Insulin Syringe with Fixed Needle	0.5mL 30G x 8mm	100
IV121127	Insulin Syringe with Fixed Needle	1mL 27G x 8mm	100
IV121129	Insulin Syringe with Fixed Needle	1mL 29G x 8mm	100
IV121130	Insulin Syringe with Fixed Needle	1mL 30G x 8mm	100
IV158029	Insulin Syringe with Retractable Safety Needle	0.5mL 29G x 8mm	100
IV158129	Insulin Syringe with Retractable Safety Needle	1mL 29G x 8mm	100
IV150127	Insulin Syringe with Retractable Safety Needle	1mL 27G x 13mm	100
IV150129	Insulin Syringe with Retractable Safety Needle	1mL 29G x 13mm	100