

Insulin Pen Needles



PRODUCT DATA SHEET

DESCRIPTION

m|devices Insulin Pen Needles consist of; a needle holder, needle tube, tip sheath, and needle seat sheath. Available in a wide range of gauges and lengths to meet all clinical requirements. This product can be used with the following brands of insulin pens: Novopen, Wanbangpen, Prandilin, HumaPen, Gansulinpen and KlikSTAR.

FEATURES AND BENEFITS

- Colour coded for easy size identification
- Gauge range from 29G - 32G
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens

PRODUCT DATA

- 32G needle is 0.23mm in 4, 6, 8mm lengths
- 31G needle is 0.25mm in 4, 5, 6, 8mm lengths
- 30G needle is 0.30mm in 6, 8, 10, 12.7mm lengths
- 29G needle is 0.33mm in 6, 8, 12, 12.7mm lengths
- ARTG Number: 272213
- Sterile EO
- Materials: PP and stainless-steel needle
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Shelf life: 5 years

RE-ORDER	DESCRIPTION	COLOUR	SIZE	CARTON/QTY
IV140431	Insulin Pen Needle 4mm	Light Blue	31G	100
IV140432	Insulin Pen Needle 4mm	Light Green	32G	100
IV140531	Insulin Pen Needle 5mm	Purple	31G	100
IV140629	Insulin Pen Needle 6mm	Pale Pink	29G	100
IV140630	Insulin Pen Needle 6mm	Brown	30G	100
IV140631	Insulin Pen Needle 6mm	Dark Blue	31G	100
IV140829	Insulin Pen Needle 8mm	Maroon	29G	100
IV140830	Insulin Pen Needle 8mm	Yellow	30G	100
IV140831	Insulin Pen Needle 8mm	Orange	31G	100
IV141030	Insulin Pen Needle 10mm	Lime Green	30G	100
IV141229	Insulin Pen Needle 12mm	Red	29G	100
IV141329	Insulin Pen Needle 10mm	Hot Pink	29G	100
IV141330	Insulin Pen Needle 12mm	Cream	30G	100