

# Insulin Pen Needles



## PRODUCT DATA SHEET

### DESCRIPTION

m|devices Insulin Pen Needles consist of; a needle holder, needle tube, tip sheath, and needle seat sheath. Available in a wide range of gauges and lengths to meet all clinical requirements. This product can be used with the following brands of insulin pens: Novopen, Wanbangpen, Prandilin, HumaPen, Gansulinpen and KlikSTAR.

### FEATURES AND BENEFITS

- Colour coded for easy size identification
- Gauge range from 29G-32G
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens

### PRODUCT DATA

- 32G needle is 0.23mm in 4mm
- 31G needle is 0.25mm in 4, 5, 6, 8mm lengths
- 30G needle is 0.30mm in 6, 8,10 lengths
- 29G needle is 0.33mm in 6, 8,12, lengths
- Materials: PP and stainless-steel needle
- Latex free
- Sterile EO
- DEHP free
- ARTG Number: 272213
- Shelf life: 5 years

RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV140431	Insulin Pen Needle	31G x 4mm	100
IV140432	Insulin Pen Needle	32G x 4mm	100
IV140531	Insulin Pen Needle	31G x 5mm	100
IV140629	Insulin Pen Needle	29G x 6mm	100
IV140630	Insulin Pen Needle	30G x 6mm	100
IV140631	Insulin Pen Needle	31G x 6mm	100
IV140829	Insulin Pen Needle	29G x 8mm	100
IV140830	Insulin Pen Needle	30G x 8mm	100
IV140831	Insulin Pen Needle	31G x 8mm	100
IV141030	Insulin Pen Needle	30G x 10mm	100
IV141229	Insulin Pen Needle	29G x 12mm	100