

Insulin Pen Needles

PRODUCT DATA SHEET



DESCRIPTION

m|devices Insulin Pen Needles consist of; a needle holder, needle tube, tip sheath, and needle seat sheath. Available in a wide range of gauges and lengths to meet all clinical requirements. This product can be used with the following brands of insulin pens: Novopen, Wanbangpen, Prandilin, HumaPen, Gansulinpen and KlikSTAR.

FEATURES AND BENEFITS

- Colour coded for easy size identification
- Gauge range from 29G - 32G
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens
- Compatible with various popular brands of insulin pens

PRODUCT DATA

- 32G needle is 0.23mm in 4, 6, 8mm lengths
- 31G needle is 0.25mm in 4, 5, 6, 8mm lengths
- 30G needle is 0.30mm in 6, 8, 10, 12.7mm lengths
- 29G needle is 0.33mm in 6, 8, 12, 12.7mm lengths
- ARTG Number: 272213
- Sterile EO
- Materials: PP and stainless-steel needle
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Shelf life: 5 years

| RE-ORDER | DESCRIPTION | COLOUR | SIZE | CARTON/QTY |
|----------|-------------------------|------------|------|------------|
| IV140531 | Insulin Pen Needle 5mm | Purple | 31G | 100 |
| IV141030 | Insulin Pen Needle 10mm | Lime Green | 30G | 100 |