

# Luer Slip Syringe without Needle

## DESCRIPTION

m|devices Luer Slip Syringes are made from polypropylene, are free from latex and are DEHP free. Available in a broad range of sizes. Available in two different nozzle positions: concentric.

## FEATURES & BENEFITS

- Bold scale markings
- Plunger stopper with double sealing ring
- Nozzle position: concentric
- Lipid resistant
- Single use only
- Peel open pouch


## PRODUCT DATA

- Available sizes: 1mL, 3mL, 5mL
- Compatible with all standard needles
- Without needle
- Sterile EO
- Material: PP
- The product is MRI suitable
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- ARTG Number: 388375
- Shelf Life: 5 Years

## ISO STANDARDS

- ISO 7886 Sterile hypodermic syringes for single use
- ISO 80369 Small-bore connectors for liquids and gases in healthcare applications
- ISO 11135:2014 Sterilisation of health-care products - Ethylene Oxide

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Luer Slip
2	Barrel
3	Plunger

## Luer Slip Syringe without Needle

STERILE

RE-ORDER	DESCRIPTION	SIZE (mL)	NOZZLE POSITION	CARTON/QTY
IV180001	Luer Slip Syringe without Needle	1mL	Concentric	100
IV180003	Luer Slip Syringe without Needle	3mL	Concentric	100
IV180005	Luer Slip Syringe without Needle	5mL	Concentric	100