

# Luer Slip Syringes



## PRODUCT DATA SHEET

### DESCRIPTION

m|devices Luer Slip Syringes are made from polypropylene, are free from latex and are DEHP free. Available in a broad range of sizes. Available in two different nozzle positions: central and eccentric

### FEATURES AND BENEFITS

- Bold scale markings
- Plunger stopper with double sealing ring
- Nozzle positions: central and eccentric
- Lipid resistant
- Single use only
- Peel open pouch

### PRODUCT DATA

- Available sizes: 1mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL and 60mL
- Compatible with all standard needles
- Without needle
- DEHP free
- Sterile EO
- Material: PP
- Shelf Life: 5 Years
- ARTG Number: 272213
- The product is MRI suitable

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV180001	Luer Slip Syringe	1mL	100
IV180003	Luer Slip Syringe	3mL	100
IV180005	Luer Slip Syringe	5mL	100
IV180010	Luer Slip Syringe	10mL	100
IV180020	Luer Slip Syringe	20mL	100
IV180030	Luer Slip Syringe	30mL	75
IV180050	Luer Slip Syringe	50mL	40
IV180060	Luer Slip Syringe	60mL	40