

# Luer Slip Syringes - Eccentric

## DESCRIPTION

m|devices Luer Slip Syringes - Eccentric, are made from polypropylene, are free from latex and are DEHP free. Available in 20mL, 30mL, 50mL and 60mL.

## FEATURES & BENEFITS

- Bold scale markings
- Plunger stopper with double sealing ring
- Nozzle positions: eccentric
- Single use only
- Peel open pouch

## PRODUCT DATA

- Available sizes: 20mL, 30mL, 50mL and 60mL
- Compatible with all standard needles
- Without needle
- Lipid resistant
- The product is MRI suitable
- Sterile EO
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Material: PP
- ARTG Number: 272213
- Shelf Life: 3 Years

## ISO STANDARDS

- ISO 7886 Sterile hypodermic syringes for single use
- ISO 80369 Small-bore connectors for liquids and gases in healthcare applications
- ISO 11135:2014 Sterilisation of health-care products - Ethylene Oxide


## PRODUCT DIAGRAM



### Nozzle Positions



Eccentric Nozzle



ITEM	DESCRIPTION
1	Luer Slip
2	Barrel
3	Flat Plunger

# Luer Slip Syringes - Eccentric

Luer Slip Syringe				STERILE
RE-ORDER	DESCRIPTION	SIZE	NOZZLE POSITION	CARTON/QTY
<b>IV180020</b>	Luer Slip Syringe	20mL	Eccentric	100
<b>IV180030</b>	Luer Slip Syringe	30mL	Eccentric	75
<b>IV180050</b>	Luer Slip Syringe	50mL	Eccentric	40
<b>IV180060</b>	Luer Slip Syringe	60mL	Eccentric	40