

# Hypodermic Syringes- Luer Lock

## PRODUCT DATA SHEET



### DESCRIPTION

m|devices luer lock three-piece hypodermic syringes are made from polypropylene, are free from latex and are DEHP free. Available in a broad range of sizes.

A retractable safety syringe version retracts the needle into the barrel after injection to prevent re-use and needle stick injury.

### FEATURES AND BENEFITS

- Bold scale markings
- Plunger stopper with double sealing ring
- Retractable version permanently disables needle reduces risk of needle stick injury and blood splatter
- Lipid resistant
- Peel open pouch
- Sterile EO
- DEHP free

### PRODUCT DATA

- Available sizes: 1mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL and 60mL
- Compatible with all standard needles
- Without needle
- Material: PP
- ARTG Number: 272213
- Shelf Life: 5 Years
- The product is MRI suitable

### PRODUCT PACKAGING LABEL



RE-ORDER	SIZE	CARTON QTY	RE-ORDER	SIZE	CARTON QTY
IV181001	1mL	100	IV181050	50mL	40
IV181003	3mL	100	IV181060	60mL	40
IV181005	5mL	100	IV170003	3mL	100
IV181010	10mL	100	IV170005	5mL	100
IV181020	20mL	100	IV170010	10mL	100
IV181030	30mL	75	IV170020	20mL	50