

# Luer Lock Syringes



## PRODUCT DATA SHEET

### DESCRIPTION

m|devices Luer Lock Syringes are three-piece hypodermic, made from polypropylene, are free from latex and DEHP free. Available in a broad range of sizes.

A retractable safety syringe version retracts the needle into the barrel after injection to prevent re-use and needle stick injury.

### FEATURES AND BENEFITS

- Bold scale markings
- Plunger stopper with double sealing ring
- Lipid resistant
- Retractable version permanently disables needle reduces risk of needle stick injury and blood splatter
- Luer slip and luer lock compatible
- Peel open pouch

### PRODUCT DATA

- Available sizes: 1mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL and 60mL
- Compatible with all standard needles
- Without needle
- The product is MRI suitable
- Material: PP
- Sterile EO
- DEHP free
- ARTG Number: 272213
- Shelf life: 5 years

RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
IV181001	Luer Lock Syringe	1mL	100
IV181003	Luer Lock Syringe	3mL	100
IV181005	Luer Lock Syringe	5mL	100
IV181010	Luer Lock Syringe	10mL	100
IV181020	Luer Lock Syringe	20mL	100
IV181030	Luer Lock Syringe	30mL	100
IV181050	Luer Lock Syringe	50mL	100
IV181060	Luer Lock Syringe	60mL	100
IV170003	Luer Lock Syringe	3mL	100
IV170005	Luer Lock Syringe	5mL	100
IV170010	Luer Lock Syringe	10mL	100
IV170020	Luer Lock Syringe	20mL	100