

Luer Activated Valve

PRODUCT DATA SHEET



DESCRIPTION

m|devices Luer Activated Valve features a closed silicone diaphragm that checks flow until mechanically opened with a luer device. Allows for bi-directional flow of up to 400mL/min upon activation. Suitable for the administration of infusion solution, drugs and for taking blood.

FEATURES AND BENEFITS

- Compatible for use with luer lock and slip connectors
- No twist necessary for luer slip connection
- Transparent housing for clear visibility
- Needleless access site with cap
- Low priming volume
- Lipid resistant

PRODUCT DATA

- Flow rate: 400mL/min
- Number of accesses: up to 100 times
- Priming volume: 0.12mL
- Pressure rating: 45 PSI
- Lipid resistant: Yes
- Negative Pressure: 284mmHg, valve opens automatically
- Material: Polycarbonate and silicone
- Sterile
- DEHP: Yes
- Latex free
- MRI compatible: Yes
- Shelf Life: 3 Years
- ARTG Number: 272212

PRODUCT PACKAGING LABEL



RE-ORDER

DESCRIPTION

CARTON/QTY

IV200000

Luer Activated Valve

50