

LuerSite Valve



PRODUCT DATA SHEET

DESCRIPTION

m|devices LuerSite Valve is designed for use with luer lock and luer slip connectors. One small push ensures a firm connection even with luer slip syringes. A closed silicone diaphragm stops flow until mechanically opened with a luer device. The LuerSite allows for a bi-directional flow of up to 345mL/min upon actuation.

FEATURES AND BENEFITS

- Needleless valve
- Bi-directional flow upon actuation
- Low priming volume
- Transparent housing for clear visibility

PRODUCT DATA

- Flow rate: 375mL/min
- Number of accesses: Tested up to 100 times
- Priming volume: 0.14mL
- Pressure rating: 45 Psi
- Lipid resistant: Yes
- Material: PC, ABS, Silicone, PE
- DEHP: Yes
- Latex free
- MRI and CT suitable: Yes
- Shelf Life: 5 Years
- ARTG Number: 272194

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	CARTON/QTY
IV200001	LuerSite Valve	50