

# Luer Slip Syringe without Needle

## DESCRIPTION

m|devices Luer Slip Syringe without Needle has a highly transparent barrel and displays clear, high contrast graduation, which allows for effortless reading and precise dosing. The conical shaped plunger stopper with double sealing ring assists with controlled aspiration and injection and allows easy elimination of air bubbles. It also features a retaining ring to prevent accidental plunger withdrawal from barrel.

Our syringes are available in a broad range of sizes: 1mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL and 60mL. They are free from latex and DEHP. It is intended to be used to draw up measured medication/solution to be delivered to the patient and aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.

## FEATURES & BENEFITS

- Conical shaped plunger stopper to assist with eliminating air bubbles from the syringe
- Highly transparent barrel with high contrast incremental measurement markings
- Double sealing ring for controlled aspiration
- Wide variety of sizes available to meet all clinical requirements


## PRODUCT DATA

- Available syringe sizes: 1mL, 3mL, 5mL, 10mL, 20mL, 30mL, 50mL and 60mL
- Without Needle
- Concentric nozzle in all size range
- Materials: PP, Polyisoprene Rubber
- Sterile: EO
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Lipid resistant: Yes
- ARTG Number: 388375
- Shelf Life: 5 Years

## ISO STANDARDS

- ISO 7886 Sterile hypodermic syringes for single use
- ISO 80369 Small-bore connectors for liquids and gases in healthcare applications
- ISO 11135:2014 Sterilisation of health-care products - Ethylene Oxide

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Luer Slip
2	Barrel
3	Plunger

# Luer Slip Syringe without Needle

Luer Slip Syringe without Needle				STERILE
RE-ORDER	DESCRIPTION	SIZE (mL)	NOZZLE POSITION	CARTON/QTY
<b>IV480001</b>	Luer Slip Syringe without Needle	1mL	Concentric	100
<b>IV480003</b>	Luer Slip Syringe without Needle	3mL	Concentric	100
<b>IV480005</b>	Luer Slip Syringe without Needle	5mL	Concentric	100
<b>IV480010</b>	Luer Slip Syringe without Needle	10mL	Concentric	100
<b>IV480020</b>	Luer Slip Syringe without Needle	20mL	Concentric	100
<b>IV480030</b>	Luer Slip Syringe without Needle	30mL	Concentric	75
<b>IV480050</b>	Luer Slip Syringe without Needle	50mL	Concentric	40
<b>IV480060</b>	Luer Slip Syringe without Needle	60mL	Concentric	40