

# Insulin Pen Needles

## 3 Bevel and 5 Bevel

### DESCRIPTION

m|devices Insulin Pen Needles are intended for use with pen injector devices for the subcutaneous injection of Insulin. Our insulin pen needles consist of seal, needle hub, needle tube, needle shield and needle container. Available in a range of gauges and lengths to meet all clinical requirements. Two needle bevels available; 3 and 5 bevel depending on the size selected.

### FEATURES & BENEFITS

- Sterile packaging
- Lubricated needle for a smooth injection
- Gauge range: 31-34G
- Needle bevel: 3 and 5 bevel needles depending on size selected
- Needle Lengths available: 4mm, 5mm, 6mm and 8mm
- Compatible with various popular brands of insulin pens, refer to IFU for complete listing.


### PRODUCT DATA

- Compatible pens: Sanofi Aventis- ClikSTAR, SoloSTAR, JuniorSTAR Lilly -KwikPen, HumaPen Memoir, HumaPen Luxura HD, HumaPen Luxura, HumaPen Sawvio, HumaPen Ergo II, Humalog Pen Novo- Nordisk FlexTouch, NovoPen3, 4 and 5, Innolet, FlexPen, Victoza, NovoPen Echo Owen Mumford -Autopen Classic, Autopen Junior Berlin Chemie- BertiPen Areo, BerlinPen Junior
- Materials: Polypropylene (PP), Polyethylene (PE), stainless steel needle, silicon oil
- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Free from PVC and Phthalates.
- Sterile EO
- ARTG Number: 392913
- Shelf Life: 5 Years

### ISO STANDARDS

- EN ISO 7864:2016 Sterile hypodermic needles for single use —Requirements and test methods
- EN ISO 11608-2:2022 Needle-based injection systems for medical use - Requirements and test methods - Part 2: Double-ended pen needles
- EN ISO 14971:2019 Medical devices - Application of risk management to medical devices

### PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Needle
2	Needle Hub

# Insulin Pen Needles

## 3 Bevel and 5 Bevel

Insulin Pen Needles 3 Bevel			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
<b>IV493315</b>	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 5mm	100
<b>IV493316</b>	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 6mm	100
<b>IV493318</b>	Insulin Pen Needle 3 Bevel	0.25mm (31G) x 8mm	100
<b>IV493324</b>	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 4mm	100
<b>IV493325</b>	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 5mm	100
<b>IV493326</b>	Insulin Pen Needle 3 Bevel	0.23mm (32G) x 6mm	75
<b>IV493334</b>	Insulin Pen Needle 3 Bevel	0.2mm (33G) x 4mm	40
<b>IV493344</b>	Insulin Pen Needle 3 Bevel	0.18mm (34G) x 4mm	40

Insulin Pen Needles 5 Bevel			STERILE
RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
<b>IV505315</b>	Insulin Pen Needle 5 Bevel	0.25mm (31G) x 5mm	100
<b>IV505316</b>	Insulin Pen Needle 5 Bevel	0.25mm (31G) x 6mm	100
<b>IV505318</b>	Insulin Pen Needle 5 Bevel	0.25mm (31G) x 8mm	100
<b>IV505324</b>	Insulin Pen Needle 5 Bevel	0.23mm (32G) x 4mm	100
<b>IV505325</b>	Insulin Pen Needle 5 Bevel	0.23mm (32G) x 5mm	100
<b>IV505326</b>	Insulin Pen Needle 5 Bevel	0.23mm (32G) x 6mm	100