

Cervical Dilating Balloon Procedure Pack

PRODUCT DATA SHEET



DESCRIPTION

m|devices Cervical Dilating Balloon Procedure Pack has been specifically designed with clinicians to support the insertion of the m|devices Cervical Dilating Balloon. The sterile pack provides all of the necessary items used for this procedure, from the sterile field for the setup trolley through to the inclusion of the m|devices Cervical Dilating Balloon, assisting the clinician with maintaining infection control principles, eliminating collection times and significantly reducing individual packaging waste.

The m|devices Cervical Dilating Balloon features a round closed distal end which meets infection control principles by maintaining a closed system and clear markings for balloon capacity on the balloon lumen and printed batch have all been considered to assist the clinician with ease of use. High contrast position markings have been printed on the shaft to assist clinicians with the placement of the Cervical Dilating Balloon.

FEATURES AND BENEFITS

- Specifically designed "fit for purpose" procedure pack
- Inclusion of the m|devices Cervical Dilating Balloon 100% Silicone, closed distal end with high contrast position markings printed on the shaft
- Significant reduction in individual packaging waste
- Pack allows for an efficient setup time and assists with maintaining infection control principles
- Batch tracking sticker provided
- Clinician's valuable time no longer spent collecting individual items

PRODUCT DATA

Contents:

1x Sterile Field 55gsm 120x120cm: CPE
 1x Blue rectangle tray:PS
 2x 2.7g Lubricating gel sterile: Glycerol and Sodium hydroxide
 4x 20mL Syringe Luer Lock: PP
 1x 25cm Sponge Forceps (Ramplcy): Stainless steel
 2x 120mL Gallipots: Modified PP

1x m|devices Cervical Dilating Balloon 80mL 20Fr 40cm: Silicone and PVC NRV
 1x XRD Gauze Swab 10x10cm 8 ply (5): cotton with X-Ray detectable filaments
 1x 15cm Clear Plastic Forceps: ABS
 1x Securement Device 8.5x 12.5cm: Non-woven fabric, velcro, acrylic adhesive and glassine paper
 1x Pouch insert: Paper

- ARTG Number: 324486
- Shelf Life: 3 years
- Latex Free
- EO Sterilisation

RE-ORDER	DESCRIPTION	CARTON/QTY
OB010002	Cervical Dilating Balloon Procedure Pack	30