

Insufflation Filter with Tubing Set

DESCRIPTION

m|devices Insufflation Filter with Tubing Sets ensure reliable and effective gas delivery for minimally invasive laparoscopic surgery. Available in 3m length and in a variety of connectors suitable for your insufflation unit. Flat filter shape also available with luer lock connectors.

FEATURES & BENEFITS

- 99.999% filtration efficiency
- Filter designed to withstand high gas flow rates and work with high-flow insufflator
- Sterile single use
- Peel open pouch
- Choice of connectors and compatible with various insufflator connectors.
- Choice of regular and flat filter



PRODUCT DATA

- Length: 3m
- Connector types: Luer connector, Push fit, O-ring, 15/22mm
- Filter- Regular and Flat shaped
- Material; PC, PVC, PP, PPH, ABS
- Not made with DEHP
- ARTG: 456226
- Shelf life: 3 years

STANDARDS

- ISO 11135:2014 Sterilisation of Healthcare products- EO
- ISO 13485:2016 Medical devices – Quality management systems
- ISO 14971:2019 – Application of risk management to medical devices

PRODUCT DIAGRAM



ITEM	DESCRIPTION	CONNECTOR TYPES
1	Luer Connector	Male Luer Lock
2	Filter	Push Fit Connector
3	Tube	O-Ring Connector
4	Connector Type	15/22mm Connector

Insufflation Filter with Tubing Set

Insufflation Filter with Tubing Set			
RE-ORDER	DESCRIPTION	LENGTH	CARTON/QTY
SU102500	Insufflation Filter and Tubing Set with Male Luer Connector, 3m smooth	3m	25
SU102501	Insufflation Filter and Tubing Set with Push Fit Connector, 3m smooth	3m	25
SU102502	Insufflation Filter and Tubing Set with O-ring Connector, 3m smooth	3m	25
SU102503	Insufflation Filter and Tubing Set with 15/22mm Connector, 3m smooth	3m	25
FLAT FILTER		LENGTH	CARTON/QTY
SU103000	Insufflation Flat Filter and Tubing Set with Male Luer Lock Connector, 3m smooth	3m	25