


2-Way Foley Silicone Catheter

PRODUCT DATA SHEET



DESCRIPTION

Indicated for both urethral and suprapubic catheterisation, our 2-Way Foley Catheter is made from 100% medical grade translucent and flexible silicone. A symmetrical balloon design supports reliable positioning and minimises the potential for ridging.

FEATURES AND BENEFITS

- Indicated for both urethral and suprapubic catheterisation
- Suitable for paediatric, male and female patients
- Symmetrical balloon
- Radiopaque
- Colour-coded
- Long indwelling time of 8 weeks
- Three lengths available: Paediatric 33cm, Female 23cm and Standard 45cm*
- Standard length codes have printed high contrast incremental markings from 6cm to 30cm*

PRODUCT DATA

- Latex free
- DEHP free
- Sterile
- Yes
- Yes
- EO
- Materials
- Shelf life
- ARTG Number

RE-ORDER	DESCRIPTION	TUBING LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
UR010006	2-Way Foley Catheter Paediatric	33cm	1.5mL	6Fr	Pink	10
UR010000	2-Way Foley Catheter Paediatric	33cm	3mL	8Fr	Light Blue	10
UR010001	2-Way Foley Catheter Paediatric	33cm	3mL	10Fr	Black	10
UR010012	2-Way Foley Catheter Female	23cm	10mL	12Fr	White	10
UR010014	2-Way Foley Catheter Female	23cm	10mL	14Fr	Green	10
UR010016	2-Way Foley Catheter Female	23cm	10mL	16Fr	Orange	10
UR011000	2-Way Foley Catheter	40cm	10mL	12Fr	White	10
UR011001	2-Way Foley Catheter	45cm	10mL	14Fr	Green	10
UR011002	2-Way Foley Catheter	40cm	10mL	16Fr	Orange	10
UR011003	2-Way Foley Catheter	40cm	10mL	18Fr	Red	10
UR011004	2-Way Foley Catheter	40cm	10mL	20Fr	Yellow	10
UR011005	2-Way Foley Catheter	45cm	10mL	22Fr	Purple	10
UR011006	2-Way Foley Catheter	45cm	10mL	24Fr	Blue	10
UR011007	2-Way Foley Catheter	40cm	10mL	26Fr	Pink	10
UR011016	2-Way Foley Catheter	40cm	10mL	16Fr	Orange	10
UR011017	2-Way Foley Catheter	40cm	30mL	18Fr	Red	10
UR011018	2-Way Foley Catheter	40cm	30mL	20Fr	Yellow	10
UR011019	2-Way Foley Catheter	40cm	30mL	22Fr	Purple	10
UR011020	2-Way Foley Catheter	40cm	30mL	24Fr	Blue	10
UR011021	2-Way Foley Catheter	40cm	30mL	26Fr	Pink	10

* 45cm stock with incremental markings is being phased-in when current 2-Way Foley Catheter stock depletes.