

2-Way Foley Catheter Tiemann Tip

DESCRIPTION

The 2-Way Foley Catheters with Tiemann Tip are made from 100% medical grade translucent and flexible silicone. Equipped with all the features of the traditional Foley Catheter, the Tiemann Tip version is designed to support the more demanding insertion. The symmetrical balloon expands equally in all directions to ensure reliable positioning in the bladder. After deflation, the balloon minimises potential ridging.

FEATURES & BENEFITS

- Symmetrical balloon
- Long indwelling time of 8 weeks
- Radiopaque
- Translucent silicone
- Colour coded
- Peel open pack


PRODUCT DATA

- Balloon size: 10mL
- Not Made with Natural Rubber Latex
- Not made with DEHP
- Materials: Silicone
- Shelf life: 3 years
- ARTG number: 211796
- Sterile EO

ISO STANDARDS

- ENISO 14155:2003 Clinical investigation for medical devices for human subjects
- ISO 10993-10:2002 Biological evaluation of the medical devices-Part 10: Tests for irritation and sensitisation
- EN550:1994 Sterilisation of medical devices-Validation and routine control of ethylene oxide sterilisation

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Side Eyelet
2	Balloon
3	X-Ray Line
4	Shaft
5	Funnel
6	Colour Code
7	Check Valve

2-Way Foley Catheter

Tiemann Tip

2-Way Silicone Foley Catheter						STERILE
RE-ORDER	DESCRIPTION	TUBING LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
UR012000	2-Way Foley Catheter	41cm	10mL	12Fr	White	10
UR012001	2-Way Foley Catheter	41cm	10mL	14Fr	Green	10
UR012002	2-Way Foley Catheter	41cm	10mL	16Fr	Orange	10
UR012003	2-Way Foley Catheter	41cm	10mL	18Fr	Red	10
UR012004	2-Way Foley Catheter	41cm	10mL	20Fr	Yellow	10
UR012005	2-Way Foley Catheter	41cm	10mL	22Fr	Purple	10
UR012006	2-Way Foley Catheter	41cm	10mL	24Fr	Blue	10
UR012007	2-Way Foley Catheter	41cm	10mL	26Fr	Pink	10