

2-Way Foley Catheter Tiemann Tip

PRODUCT DATA SHEET



DESCRIPTION

The 2-Way Foley Catheters with Tiemann Tip are made from 100% medical grade translucent and flexible silicone. Equipped with all the features of the traditional Foley Catheter, the Tiemann Tip version is designed to support the more demanding insertion. The symmetrical balloon expands equally in all directions to ensure reliable positioning in the bladder. After deflation, the balloon minimises potential ridging.

FEATURES AND BENEFITS

- Radiopaque
- Symmetric Balloon
- Long indwelling time of 8 weeks
- Translucent silicone
- Colour coded
- Peel open pack

PRODUCT DATA

- Length: 41cm
- Balloon size: 10mL
- DEHP free
- Latex free
- Materials: Silicone
- Shelf life: 3 years
- ARTG number: 211796
- Sterile EO

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
UR012000	2-Way Foley Catheter TT	41cm	10mL	12Fr	White	10pcs
UR012001	2-Way Foley Catheter TT	41cm	10mL	14Fr	Green	10pcs
UR012002	2-Way Foley Catheter TT	41cm	10mL	16Fr	Orange	10pcs
UR012003	2-Way Foley Catheter TT	41cm	10mL	18Fr	Red	10pcs
UR012004	2-Way Foley Catheter TT	41cm	10mL	20Fr	Yellow	10pcs
UR012005	2-Way Foley Catheter TT	41cm	10mL	22Fr	Purple	10pcs
UR012006	2-Way Foley Catheter TT	41cm	10mL	24Fr	Blue	10pcs
UR012007	2-Way Foley Catheter TT	41cm	10mL	26Fr	Pink	10pcs