

# 3-Way Silicone Foley Catheter

## PRODUCT DATA SHEET



### DESCRIPTION

The 3-Way Silicone Foley Catheter is made from 100% medical grade translucent and flexible silicone. The symmetrical balloon design supports reliable positioning and minimises for potential ridging.

### FEATURES AND BENEFITS

- Urethral and SPC approved
- Symmetrical balloon
- Long indwelling time 30 days
- Radiopaque
- Colour coded
- Peel open pack

### PRODUCT DATA

- Latex free
- DEHP free
- Materials: Silicone
- Shelf life: 36 years
- ARTG number: 211796
- Sterile EO

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
UR013000	3-Way Foley Catheter	40cm	30mL	16Fr	Orange	10pcs
UR013001	3-Way Foley Catheter	43cm	30mL	18Fr	Red	10pcs
UR013002	3-Way Foley Catheter	43cm	50mL	20Fr	Yellow	10pcs
UR013003	3-Way Foley Catheter	43cm	50mL	22Fr	Purple	10pcs
UR013004	3-Way Foley Catheter	43cm	50mL	24Fr	Blue	10pcs