

# 3-Way Foley Catheter Silicone

## DESCRIPTION

The 3-Way Silicone Foley Catheter is made from 100% medical grade translucent and flexible silicone and adds continuous irrigation capabilities for use in removing debris and blood clots. The symmetrical balloon design supports reliable positioning and minimises for potential ridging.

## FEATURES & BENEFITS

- Symmetrical balloon
- Long indwelling time 29 days
- Radiopaque
- Colour coded
- Peel open pack


## PRODUCT DATA

- Not Made with Natural Rubber Latex
- DEHP NO = Not detected (Result less than the Method Detection Limit EN 14372)
- Materials: Silicone
- Shelf life: 3 years
- ARTG number: 211796
- Sterile EO

## ISO STANDARDS

ISO 14971 - Risk Management to Medical Devices  
 ISO 80369 - Small bore Connectors  
 ISO 20696 - Sterile Urethral Catheters for single use

## PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Side Eyelet
2	Balloon
3	X-Ray Line
4	Shaft
5	Funnel
6	Colour Code
7	Check Valve

# 3-Way Foley Catheter Silicone

3-Way Silicone Foley Catheter						STERILE
RE-ORDER	DESCRIPTION	TUBING LENGTH	BALLOON CAPACITY	GAUGE	COLOUR	CARTON/QTY
<b>UR013000</b>	3-Way Foley Catheter	40cm	30mL	16Fr	Orange	10
<b>UR013001</b>	3-Way Foley Catheter	43cm	30mL	18Fr	Red	10
<b>UR013002</b>	3-Way Foley Catheter	43cm	50mL	20Fr	Yellow	10
<b>UR013003</b>	3-Way Foley Catheter	43cm	50mL	22Fr	Purple	10
<b>UR013004</b>	3-Way Foley Catheter	43cm	60mL	24Fr	Blue	10