

Cystoscopy Irrigation Set

For Cystoscopy/TUR Irrigation

DESCRIPTION

The m|devices Cystoscopy TUR Irrigation Set is designed for Cystoscopy or Bladder Irrigation after Transurethral resection (TUR) or bladder tumour. It is made from medical grade non-toxic PVC and utilises a soft flexible drip chamber and delivers high flow rates.

FEATURES & BENEFITS

- Double spike configuration
- Soft chamber without filter
- Stepped connector with silicone tube end
- Two click clamps


PRODUCT DATA

- Tubing inner diameter: 4.8mm
- Tubing outer diameter: 6.8mm
- Priming volume: 72mL
- Sterile EO
- Material: PVC
- ARTG Number: 314960
- Not Made with Natural Rubber Latex
- DEHP YES = Product contains Phthalates such as the PVC Plasticizer Di-(2-ethylhexyl)-phthalate (DEHP).
- Shelf life: 5 Years

STANDARDS

- EN ISO 10993-5:2009 Biological evaluation of medical devices - Tests for in vitro cytotoxicity.
- EN ISO 10993-7:2008 Biological evaluation of medical devices - Ethylene oxide sterilisation residuals
- EN ISO 14971:2012 Application of Risk Management.

PRODUCT DIAGRAM



ITEM	DESCRIPTION	6	Tube
1	Cover	7	Roller Body
2	Spike Irrigation	8	Roller wheel
3	Clamp	9	Stepped connector
4	Y-connector	10	Tube
5	Soft Chamber	11	Tube

Cystoscopy TUR Irrigation Set

RE-ORDER	DESCRIPTION	TOTAL LENGTH	CARTON/QTY
UR020001	Cystoscopy TUR Irrigation Set	285cm	50