


Nelaton Catheters Standard



PRODUCT DATA SHEET

DESCRIPTION

Our Nelaton Catheters are intended for intermittent bladder catheterisation. They are made from medical grade thermo-sensitive plastic and offer a perfectly finished distal end for smooth insertion.

FEATURES AND BENEFITS

- Thermo-sensitive PVC
- Heat polished eyelets
- Smooth round tip
- Colour coded connectors
- Available for Paediatric, Female and Male patients

PRODUCT DATA

- Latex free
- DEHP free
- Materials: PVC
- Shelf life: 5 Years
- ARTG Number: 210992

RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
UR040000	Standard Nelaton Catheter Paediatric	30cm	6Fr	Light Green	50
UR040001	Standard Nelaton Catheter Paediatric	30cm	8Fr	Light Blue	50
UR040002	Standard Nelaton Catheter Paediatric	30cm	10Fr	Black	50
UR041000	Standard Nelaton Catheter Female	20cm	8Fr	Light Blue	50
UR041001	Standard Nelaton Catheter Female	20cm	10Fr	Black	50
UR041002	Standard Nelaton Catheter Female	20cm	12Fr	White	50
UR041003	Standard Nelaton Catheter Female	20cm	14Fr	Green	50
UR041004	Standard Nelaton Catheter Female	20cm	16Fr	Orange	50
UR041005	Standard Nelaton Catheter Female	20cm	18Fr	Red	50
UR041006	Standard Nelaton Catheter Female	20cm	20Fr	Yellow	50
UR042000	Standard Nelaton Catheter Male	40cm	8Fr	Light Blue	50
UR042001	Standard Nelaton Catheter Male	40cm	10Fr	Black	50
UR042003	Standard Nelaton Catheter Male	40cm	12Fr	White	50
UR042003	Standard Nelaton Catheter Male	40cm	14Fr	Green	50
UR042004	Standard Nelaton Catheter Male	40cm	16Fr	Orange	50
UR042005	Standard Nelaton Catheter Male	40cm	18Fr	Red	50
UR042006	Standard Nelaton Catheter Male	40cm	20Fr	Yellow	50