


Hydrophilic Coated Nelaton Catheters



PRODUCT DATA SHEET

DESCRIPTION

Our Hydrophilic Coated Nelaton Catheters are made from non-toxic DEHP free PVC. Their smooth water-soluble polymer coating combined with a heat polished distal tip supports a smooth insertion without the need of lubricants.

FEATURES AND BENEFITS

- Thermo-sensitive PVC
- Heat polished eyelets
- Smooth round tip
- Colour coded connectors
- Available for Paediatric, Female and Male patients
- Material: PVC
- DEHP: Yes
- Latex: No
- ARTG Number: 210992
- Shelf Life: 5 Years

PRODUCT DATA

RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON QTY
UR043001	Hydrophilic Coated Nelaton Catheter Paediatric	30cm	8Fr	Light Blue	25
UR043002	Hydrophilic Coated Nelaton Catheter Paediatric	30cm	10Fr	Black	25
UR044000	Hydrophilic Coated Nelaton Catheter Female	18cm	8Fr	Light Blue	25
UR044001	Hydrophilic Coated Nelaton Catheter Female	18cm	10Fr	Black	25
UR044002	Hydrophilic Coated Nelaton Catheter Female	18cm	12Fr	White	25
UR044003	Hydrophilic Coated Nelaton Catheter Female	18cm	14Fr	Green	25
UR044004	Hydrophilic Coated Nelaton Catheter Female	18cm	16Fr	Orange	25
UR044005	Hydrophilic Coated Nelaton Catheter Female	18cm	18Fr	Red	25
UR045000	Hydrophilic Coated Nelaton Catheter Male	40cm	8Fr	Light Blue	25
UR045001	Hydrophilic Coated Nelaton Catheter Male	40cm	10Fr	Black	25
UR045002	Hydrophilic Coated Nelaton Catheter Male	40cm	12Fr	White	25
UR045003	Hydrophilic Coated Nelaton Catheter Male	40cm	14Fr	Green	25
UR045004	Hydrophilic Coated Nelaton Catheter Male	40cm	16Fr	Orange	25
UR045005	Hydrophilic Coated Nelaton Catheter Male	40cm	18Fr	Red	25