

# Hydrophilic Nelaton Catheter with Water Sachet -Paediatric

## PRODUCT DATA SHEET



### DESCRIPTION

m|devices Hydrophilic Nelaton Catheter with Water Sachet in paediatric length, also contains a blue coloured holding sleeve for easy visualisation for an aseptic insertion. This all-in-one product which offers convenience and ease of use, without the need for additional products.

The hydrophilic coating, heat polished closed distal end and bilateral side eyelets support a smooth and painless insertion and efficient drainage of urine.

### FEATURES AND BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Smooth round tip
- Blue coloured holding sleeve to maintain sterility during insertion
- Colour coded connectors
- Latex free
- Heat polished eyelets

### PRODUCT DATA

- Paediatric length: 30cm
- Sizes 6Fr-10Fr
- DEHP: Yes
- Shelf Life: 3 Years
- Material: Nelaton: PVC
- Sachet: Purified Water
- Holding Sleeve: PE
- ARTG Number: 202850

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
UR046000	Hydrophilic Nelaton with Water Sachet - Paediatric	30cm	6Fr	Light Green	25
UR046001	Hydrophilic Nelaton with Water Sachet - Paediatric	30cm	8Fr	Light Blue	25
UR046002	Hydrophilic Nelaton with Water Sachet - Paediatric	30cm	10Fr	Black	25