

Hydrophilic Nelaton with Water Sachet Paediatric

PRODUCT DATA SHEET



DESCRIPTION

m|devices Hydrophilic Nelaton Catheter with Water Sachet provides an all-in-one product which offers convenience and ease of use, without the need for additional products. The hydrophilic coating, heat polished closed distal end and bilateral side eyelets support a smooth and painless insertion and efficient drainage of urine.

FEATURES AND BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Blue coloured holding sleeve to maintain sterility during insertion
- Colour coded connectors
- Heat polished eyelets
- Smooth round tip

PRODUCT DATA

- Paediatric length: 30cm
- Sizes 6Fr-10Fr
- DEHP free
- Latex free
- Material: Nelaton: PVC
- Sachet: Purified water
- Holding sleeve: PE
- Shelf life: 3 years
- ARTG number: 210992

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
UR046000	Hydrophilic Nelaton with Water Sachet - Paediatric	30cm	6Fr	Light Green	25pcs
UR046001	Hydrophilic Nelaton with Water Sachet - Paediatric	30cm	8Fr	Light Blue	25pcs
UR046002	Hydrophilic Nelaton with Water Sachet - Paediatric	30cm	10Fr	Black	25pcs