


Hydrophilic Nelaton with Water Sachet Female

PRODUCT DATA SHEET



DESCRIPTION

m|devices Hydrophilic Nelaton Catheter with Water Sachet provides an all-in-one product which offers convenience and ease of use, without the need for additional products. The hydrophilic coating, heat polished closed distal end and bilateral side eyelets support a smooth and painless insertion and efficient drainage of urine.

FEATURES AND BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Smooth round tip
- Colour coded connectors
- Heat polished eyelets
- Blue coloured holding sleeve to maintain sterility during insertion

PRODUCT DATA

- Female length: 18cm
- Sizes 8Fr- 18Fr
- Material: Nelaton: PVC
- Holding Sleeve: PE
- Sachet: Purified Water
- DEHP free
- Latex free
- Shelf Life: 3 Years
- ARTG Number: 210992

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
UR047000	Hydrophilic Nelaton with Water Sachet - Female	18cm	8Fr	Light Blue	25pcs
UR047001	Hydrophilic Nelaton with Water Sachet - Female	18cm	10Fr	Black	25pcs
UR047002	Hydrophilic Nelaton with Water Sachet - Female	18cm	12Fr	White	25pcs
UR047003	Hydrophilic Nelaton with Water Sachet - Female	18cm	14Fr	Green	25pcs
UR047004	Hydrophilic Nelaton with Water Sachet - Female	18cm	16Fr	Orange	25pcs
UR047005	Hydrophilic Nelaton with Water Sachet - Female	18cm	18Fr	Red	25pcs