

Hydrophilic Nelaton with Water Sachet Male

PRODUCT DATA SHEET



DESCRIPTION

m|devices Nelaton Catheter with Water Sachet in a male length, also contains a blue coloured holding sleeve for easy visualisation for an aseptic insertion. This all in one product offers convenience and ease of use without the need for additional products.

The hydrophilic coating, heat polished closed distal end and bilateral side eyelets support a smooth and painless insertion and efficient drainage of urine.

FEATURES AND BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Smooth round tip
- Colour coded connectors
- Heat polished eyelets
- Blue coloured holding sleeve to maintain sterility during insertion

PRODUCT DATA

- Male length: 40cm
- Sizes 8Fr-18Fr
- Material: Nelaton: PVC
- Sachet: Purified Water
- Holding Sleeve: PE
- DEHP free
- Latex free
- Shelf Life: 3 Years
- ARTG Number: 210992

PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	LENGTH	GAUGE	COLOUR	CARTON/QTY
UR048000	Hydrophilic Nelaton with Water Sachet - Male	40cm	8Fr	Light Blue	25pcs
UR048001	Hydrophilic Nelaton with Water Sachet - Male	40cm	10Fr	Black	25pcs
UR048002	Hydrophilic Nelaton with Water Sachet - Male	40cm	12Fr	White	25pcs
UR048003	Hydrophilic Nelaton with Water Sachet - Male	40cm	14Fr	Green	25pcs
UR048004	Hydrophilic Nelaton with Water Sachet - Male	40cm	16Fr	Orange	25pcs
UR048005	Hydrophilic Nelaton with Water Sachet - Male	40cm	18Fr	Red	25pcs