

Hydrophilic Nelaton with Water Sachet Coude Tip

PRODUCT DATA SHEET



DESCRIPTION

m|devices Hydrophilic Nelaton Catheter with Water Sachet Coude Tip contains a blue coloured holding sleeve for easy visualisation for an aseptic insertion. This all in one product offers convenience and ease of use without the need for additional products.

The hydrophilic coating in conjunction with the Coude angled closed distal end and bilateral side eyelets, support a smooth and painless insertion and efficient drainage of urine.

FEATURES AND BENEFITS

- Hydrophilic coated
- All-in-one product with water sachet
- Smooth round tip
- Coude tip to support difficult insertions
- Colour coded connectors
- Heat polished eyelets
- Blue coloured holding sleeve

PRODUCT DATA

- Length: 40cm
- Sizes 8Fr- 18Fr
- Material: Nelaton: PVC
- Sachet: Purified Water
- Holding Sleeve: PE
- DEHP free
- Latex free
- Shelf Life: 3 Years
- ARTG Number: 210992

PRODUCT PACKAGING LABEL



| RE-ORDER | DESCRIPTION | LENGTH | GAUGE | COLOUR | CARTON/QTY |
|----------|---|--------|-------|--------|------------|
| UR049001 | Hydrophilic Nelaton with Water Sachet - Coude Tip | 40cm | 10Fr | Black | 25pcs |
| UR049002 | Hydrophilic Nelaton with Water Sachet - Coude Tip | 40cm | 12Fr | White | 25pcs |
| UR049003 | Hydrophilic Nelaton with Water Sachet - Coude Tip | 40cm | 14Fr | Green | 25pcs |
| UR049004 | Hydrophilic Nelaton with Water Sachet - Coude Tip | 40cm | 16Fr | Orange | 25pcs |
| UR049005 | Hydrophilic Nelaton with Water Sachet - Coude Tip | 40cm | 18Fr | Red | 25pcs |