

Universal Securement Device

DESCRIPTION

The Universal Securement Device provides securement of catheters, drains and other tubing in two simple steps. Adhesive strip ensures the catheter does not slide once in position and the soft material tabs cross over each other to secure the tubing in place. Available in two sizes; standard and large, to accommodate a variety of diameters of tubing, catheters, drains and other devices.

FEATURES & BENEFITS

- Two stabilisation steps ensures maximum securement to prevent sliding and dislodgement
- Low irritant, soft adhesive dressing
- Soft material tabs can be re-applied numerous times to reposition
- Easy to apply and remove
- Upon removal skin remains free from residue and does not require any additional products to assist with removal
- Two available sizes accommodates various tubing diameters and provides versatility


PRODUCT DATA

- Standard size: 6.5cm x 10.5cm
- Large size: 8.5cm x 12.5cm
- MRI compatible: Yes
- Materials: Nylon, non-woven fabric, acrylic adhesive, glassine paper
- Not Made with Natural Rubber Latex
- ARTG Number: 299241
- Shelf life: 3 Years

ISO STANDARDS

- ISO 11135:2014 Sterilisation of health care products -- Ethylene oxide -- Requirements for development, validation and routine control of a sterilization process for medical devices.
- ISO 11138-1:2017 Sterilisation of health care products —Biological indicators

PRODUCT DIAGRAM



ITEM	DESCRIPTION
1	Tabs
2	Adhesive Strip
3	Fabric Fastener

Universal Securement Device

STERILE

RE-ORDER	DESCRIPTION	SIZE	CARTON/QTY
UR060000	Universal Securement Device	6.5cm x 10.5cm	50
UR060001	Universal Securement Device	8.5cm x 12.5cm	50