

# Drainage Bags

## PRODUCT DATA SHEET



### DESCRIPTION

The Drainage Bags are used to collect and contain urine via an indwelling catheter or leg bag. Our range of urine bags are available in sterile and non-sterile versions with and without a non-return valve. Variations in tube length ensure the bag is optimised for different settings.

### FEATURES AND BENEFITS

- Durable material
- Clear volume graduation
- With or without non-return valve
- With T-Tap
- With or without needleless sampling port
- Varying tube lengths to support patient's needs
- Latex free

### PRODUCT DATA

- Material: PVC, PP, PE
- DEHP: Yes
- Sterile: EO
- ARTG Number: 267875
- Shelf Life: 5 Years

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	TUBING ID	TUBING OD	CARTON QTY
UR062000	Urine bag 2000mL T-TAP NRV	90cm	5mm	6.3mm	250
UR062001	Urine bag 2000mL T-TAP	90cm	5mm	6.3mm	250
UR062002	Urine bag 2000mL T-TAP NRV Frosted Tubing Sample Port	100cm	6.5mm	9.3mm	100
UR062003	Urine bag 2000mL T-TAP NRV Clear Tubing Sample Port	120cm	7.4mm	9.8mm	100