

# Urine Bags

## PRODUCT DATA SHEET



### DESCRIPTION

The Urine Bags are used to collect and contain urine via an indwelling catheter or leg bag. Our range of Urine Bags are available in sterile and non-sterile versions with and without a non-return valve. Variations in tube length ensure the bag is optimised for different settings.

### FEATURES AND BENEFITS

- Durable material
- Clear volume graduation
- With or without non-return valve
- With T-Tap
- With or without needleless sampling port
- Varying tube lengths to support patient's needs

### PRODUCT DATA

- Material: PVC, PP, PE
- DEHP: Yes
- Latex free
- Sterile EO
- ARTG Number: 267875
- Shelf life: 5 Years

### PRODUCT PACKAGING LABEL



RE-ORDER	DESCRIPTION	TUBING LENGTH	BAG CAPACITY	CARTON/QTY
UR062000	Urine Bag T-Tap NRV ID 5mm OD 6.3mm	90cm	2000mL	250pcs
UR062001	Urine Bag T-Tap ID 5mm OD 6.3mm	90cm	200mL	250pcs
UR062002	Urine Bag T-Tap NRV, Sample Port, Frosted Tube, ID 6.5mm OD 9.3mm	100cm	2000mL	100pcs
UR062003	Urine Bag T-Tap NRV, Sample Port, Clear Tube, ID 7.4mm OD 9.8mm	120cm	2000mL	100pcs
UR062007	Urine bag 2000mL T-TAP NRV	110cm	2000mL	250pcs
UR062008	Urine Bag T-Tap	110cm	2000mL	250pcs